

**Specifications :-**

- To Smoothen Amalgam After Condensing, Used To Create Occlusal Anatomy, Burnish Amalgam
- Acorn Shaped Tip, Metal, Double Sided.
- Standard Stainless Steel
- Handle
- Autoclaveable.

*Note: Specifications are subject to change, Photos shown above are Indicative, Actual Product can Vary.*



Export Sales: +91-9829132777  
India Sales: +91-9588842361



IT-2013, Ramchandrapura Industrial Area,  
Sitapura Extension, Jaipur-302022, India.



info@tesca.in  
www.tescaglobal.com