



## **Applications:**

Used to determine initial and final set times of portland cement, masonary cement, hydraulic hydrated lime and certain mortars.

## **Details:**

- The adjustable support has horizontal arms to guide two weighted needles with flat-end cylindrical stainless steel tips.
- The support assembly is adjustable. Base has ample flat shelf for positioning specimens.
- The two steel needle weights are calibrated to meet specifications.
- The needle points are manufactured from stainless steel.
- The initial setting needle is 2.12 mm in diameter and weighs 113 grams.
- The final setting needle is 1.06 mm in diameter and weighs 453.6 grams.
- Supplied with 02 Nos glass plates Dimension  $100 \pm 5$  mm square
- Also supplied complete with 02 Nos Cement paste pat Mould (Base diameter 76  $\pm$  13 mm & Top diameter  $50 \pm 13 \, mm$ )

## Features:

- Meets several ASTM and AASHTO test methods for set-time testing
- · Includes initial and final-set needles
- Large base to position specimens

Overall Dimensions: 230 X 100 X 260 (L X W X D) mm

Note: Specifications are subject to change.

Near Bombay Hospital, Vidhani Circle, Jaipur-302022, Rajasthan, India,

Tel: +91-9829132777; Email: info@tesca.in, tesca.technologies@gmail.com

<sup>™</sup> Website: www.tescaglobal.com

