

**Details**

- Loading pad made of perspex for 50mm dia for specimens.
- Pair of plain discs made of perspex for 50mm dia X 6mm thick.
- Pair of porous stone for 50mm dia X 6mm thick.
- Split sand former for 50mm dia specimen.
- Sheath stretcher for 50mm dia specimen.
- One dozen rubber sheaths for 50mm dia samples.
- For synthetic rubber 'O' rings for 50mm dia specimen.
- One Plastic drainage tube.
- Overall Dimensions : 240 X 240 X 340 (L X W X D) mm

Note: Specifications are subject to change.